



Edition: 11.22
Valid until: 12.23
Rev. Jan 22



XENON² SA

Deliver to:

Name:

Street:

Town:

E-Mail:

Tel:

Postal code:

Fax:



Deliver to:

Name:

Street:

Town:

E-Mail:

Tel:

Postal code:

Fax:

Order date:

Purchase order:

Number required:

☐ Order ☐ Calculation ☐ Built 4 me (longer leadtime)

Art. No. DE: 77 800 000

Marked for:

Sunrise order no:

● = Standard

○ = No cost option

☐ = Cost option

Max. user weight:

125 kg

USER WEIGHT

kg

Mandatory information to ensure correct configuration

Base Price

Xenon² SA

● = Standard ○ = No cost option ☐ = Cost option

CLIENT MEASUREMENTS (all measurements without cushion)

SEAT WIDTH SELECTION - measured from outside frame tube to outside frame tube (mm)

320	340	360	380	400	420	440	460	480	500
<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>

SEAT DEPTH SELECTION - measured in mm

340	360	380	400	420	440	460	480	500
<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
380 mm frame length			each frame in a separate length (20 mm increments)					

FRONT SEAT HEIGHT SELECTION: in dependance of front castors, Carbotecture forks and frame versions

		Front seat height in mm															
		Fork (length)	430	440	450	460	470	480	490	500	510	520	530	540	550	560	
Front Castor	3"	111 mm	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>												
		123 mm (1)		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>											
		138 mm (1)				<input type="radio"/>	<input type="radio"/>	<input type="radio"/>									
		174 mm							<input type="radio"/>	<input type="radio"/>	<input type="radio"/>						
		Frog Leg (1)					<input type="checkbox"/>	<input type="checkbox"/>									
		One Arm (1)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>												
	4"	111 mm	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>											
		123 mm (1)		<input type="checkbox"/>	<input type="checkbox"/>		<input type="checkbox"/>	<input type="checkbox"/>									
		138 mm (1)				<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>								
		174 mm								<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>				
		Frog Leg (1)					<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>								
		One Arm (1)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>											
	5"	111 mm			<input type="radio"/>		<input type="radio"/>	<input type="radio"/>									
		123 mm (1)					<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>								
		138 mm (1)					<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>							
		174 mm							<input type="radio"/>	<input type="radio"/>		<input type="radio"/>	<input type="radio"/>	<input type="radio"/>			
		Frog Leg (1)						<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>							
		One Arm (1)			<input type="checkbox"/>		<input type="checkbox"/>	<input type="checkbox"/>									
		123 mm (1)						<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>							
		138 mm (1)						<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>						
174 mm									<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>			
Frog Leg (1)									<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>						
7"		138 mm (1)							<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>					
	174 mm										<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>			

REAR SEAT HEIGHT SELECTION: in dependance of rear wheel size

		Rear seat height in mm													
Rear Wheels		370	380	390	400	410	420	430	440	450	460	470	480	490	500
24"		<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
25"		<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>

The rear seat height can only be between 0 mm and 110 mm below the selected front seat height.

(1) approved for use with powered add-ons (powered wheels or push-devices).

Please be advised that powered add-ons can not be approved in combination with 111mm-forks and 174mm-forks. In general, always verify whether a combination agreement is in place before fitting powered add-ons.

(2) only with lightweight adapter

For more information
visit our webpage

www.SunriseMedical.eu

© Sunrise Medical 2023

Please follow this measuring method:

The following method is recommended on the basis of extensive practical experience. Deviations are possible in individual instances.

Seat width (SW):

- a) For the user: Allow space between sideguard and hip on each side of the chair
(approx. 10 - 20 mm altogether)
- b) On the wheelchair: From outside frame tube to outside frame tube.

Seat depth (SD):

- a) For the user: When sitting upright: leave space for 2 - 3 fingers between
front edge of the seat upholstery and back of knee.
(approx. 30 - 60 mm). Care must be taken with adj. upholstery.
- b) On the wheelchair: From back tube to front of seat upholstery.

Front seat height:

- a) For the user: When sitting with 90° kneeangle:
From back of knee to bottom of shoe sole.
- b) On the wheelchair: From front of seat tube to the ground.
(min. 30 mm longer than user measurement)
Remember to include the cushion height.

Rear seat height:

- b) On the wheelchair: From back of seat tube to the ground.
Note: Max. seat height difference rear to front is 110 mm

Back height:

- a) For the user: From the seat, including cushion, to be no more than 10 - 20 mm below the lower shoulder blade.
- b) On the wheelchair: From back of seat tube to upper edge of backrest upholstery.
Note: The user's ability to sit upright has a major influence on the
back height and it can therefore differ considerably

Back angle:

Important: The back height depends on the degree of disability.
0° is the angle 90° to the ground. Negative angles are <90°
and positive angles are >90°.

Technical Data Overview:

Max. overall height:	1010 mm		
Max. overall length:	1150 mm		
Max. turning circle:	720 mm		
Max. safe slope:	10°		
Wheelchair weight:	from 10,3 kg		
Dismantled weight:	from 7,8 kg		
Weight of heaviest part:	2,1 kg		
Max. overall width:	Camber	24" rear wheel	25" rear wheel
	For 0° camber:	SW + 160 mm	SW + 160 mm
	For 2° camber:	SW + 210 mm	SW + 220 mm
	For 4° camber:	SW + 260 mm	SW + 270 mm

Product is **crash tested according to ISO 7176-19** and therefore approved for transportation in a vehicle.
in combination with XES090100 (Tie down brackets for transportation)
For transport without a user the wheelchair must be safely clamped in the vehicle.

All changes in the configuration might change the seat heights.
Manufacturing tolerances of +/- 10 mm

Xenon² SA

● = Standard ○ = No cost option □ = Cost option

Price in €

SEAT

Seat sling

XES020002	●	Seat sling, nylon, tension adjustable, 1 utility bag	
XES020003	□	Seat sling, nylon, tension adjustable, 2 utility bags (only from SW 360 mm)	
XES020004	□	Seat sling, nylon, tension adjustable, 1 utility bag, 1 catheter bag (only from SW 360 mm)	
XES020000	○	Without seat sling	

Seat cushion, cover with zipper, black

Foam (medium-soft)	Thickness:	XFF020101	□ 30 mm	XFF020102	□ 50 mm	
Latex (very soft)	Thickness:	XFF020103	□ 50 mm			

Seat accessories

XES090018	□	Lap belt with buckle, metal	
XES020050	□	Additional padding for seat upholstery, removable	

JAY Cushions - for complete product range see Jay order forms

JAY020002	□	JAY Basic	
JAY020003	□	JAY Soft Combi P	
JAY020001	□	JAY Easy Visco, Curved Base	
JAY020004	□	JAY Lite	
JAY020007	□	JAY Easy Fluid, Curved Base	

BACKREST

Back height range

Choose back height factory setting (mm).

Height adjustments of in total 50 mm possible.

250 - 300	300 - 350	350 - 400	400 - 450	450 - 475
○ 250	○ 300	○ 350	○ 400	○ 450
○ 275	○ 325	○ 375	○ 425	○ 475
○ 300	○ 350	○ 400	○ 450	

Backrest height, carbon fiber tubes (1)

Fixed height, non height adjustable

□ 250	□ 275	□ 300	□ 325
□ 350	□ 375	□ 400	□ 425

XES030010	○	Backrest, fixed (90°)	
XES030029	○	Desk sideguard preparation for later replacement (2)	
XES030031	□	Backrest, fold-down (to the front) (4) (only available from back height 325 mm upwards)	
XES030012	□	Backrest, angle-adjustable (75° to 103°)	
XES030029	○	Desk sideguard preparation for later replacement (3)	
XES030031	□	Backrest, fold-down (to the front) (4) (only available from back height 325 mm upwards)	

Pushhandles

XES030200	○	Without push handles	
XES030201	○	Push handles long	
XES030203	□	Push handles, fold-down (only available from back height 325 mm upwards)	
XES030205	□	Push handles, height-adjustable for Jay backs (only available from back height 325 mm upwards)	
XES030204	□	Push handles, height-adjustable (only available from back height 300 mm upwards)	

Backrest upholstery

XES030316	○	EXO Backrest sling, padded, adjustable, 1 accessory bag, black (EXO Evo version) Colour, binding tape: ● Black ○ Orange ○ Silver ○ Red ○ Blue	
XES030317	□	EXO PRO Backrest, padded, adj., breathable spacer fabric, 1 accessory bag, black (EXO Evo version) Colour, binding tape: ● Black ○ Orange ○ Silver ○ Red ○ Blue	
XES030300	○	Without backrest upholstery	


Accessories for backrest

XES030402	□	Stabilising bar, automatic folding (only available from SW 380mm; not with fixed height backrest)	
XES030401	□	Stabilising bar, swing-away (only available from SW380) □ Headrest, height-, depth- and angle-adjustable XES030411 ○ medium (200x110mm) XES030412 ○ tall (230x110mm)	

- (1) not with Jay J3 backs, not with Desk sideguard
 (2) necessary for retrofitting Desk sideguard or Tubular armrest with Backrest, fixed 90°
 (3) necessary for retrofitting Desk sideguard
 (4) not with height adj. push handles, not with JAY backs; in combination with push handles fold-down
 only from backrest height 400mm upwards; please see Folding Matrix on next side

For more information
 visit our webpage
www.SunriseMedical.eu

© Sunrise Medical 2023



XES 020002 Seat sling with utility bag 1



XES 020003 Seat sling with utility bags 2



XES 030316 EXO backrest sling (EXO Evo version)



XES 030317 EXO PRO backrest sling (EXO Evo version)



XES 030201 Push handles long



XES 030203 Push handles, fold-down



XES 030204 Push handles, height-adjustable



XES 090018 Lap belt with buckle, metal

Folding-Matrix for XES 030031 Backrest, fold-down to the front.



Seat width

Back height

	325	350	375	400	425	450	475
320	✓	✓	✓	✓	✓	✓	✓
340	✓	✓	✓	✓	✓	✓	✓
360	✓	✓	✓	✓	✓	✓	✓
380		✓	✓	✓	✓	✓	✓
400			✓	✓	✓	✓	✓
420				✓	✓	✓	✓
440				✓	✓	✓	✓
460					✓	✓	✓
480						✓	✓
500							✓



=folds 90° to the front on folded chair



=can't be folded full 90° to the front on folded chair in combination with fold-down pushhandles



=can't be folded full 90° to the front on folded chair



=will fold less than 45° to the front on folded chair, folds 90° on unfolded chair

Xenon² SA

● = Standard ○ = No cost option □ = Cost option

Price in €

Jay 3 Back - Shallow Contour

Mid Thoracic:	<input type="checkbox"/> JAY030010	Medium - 340 mm	<input type="checkbox"/> JAY030033	Tall - 380 mm	
Upper Thoracic:	<input type="checkbox"/> JAY030012	Short - 420 mm	<input type="checkbox"/> JAY030019	Medium - 460 mm	<input type="checkbox"/> JAY030025 Tall - 500 mm
Shoulder Height:	<input type="checkbox"/> JAY030015	Short - 530 mm	<input type="checkbox"/> JAY030022	Medium - 570 mm	<input type="checkbox"/> JAY030028 Tall - 610 mm

Jay 3 Back - Mid Contour

Mid Thoracic:	<input type="checkbox"/> JAY030011	Medium - 340 mm	<input type="checkbox"/> JAY030018	Tall - 380 mm	
Upper Thoracic:	<input type="checkbox"/> JAY030013	Short - 420 mm	<input type="checkbox"/> JAY030020	Medium - 460 mm	<input type="checkbox"/> JAY030026 Tall - 500 mm
Shoulder Height:	<input type="checkbox"/> JAY030016	Short - 530 mm	<input type="checkbox"/> JAY030023	Medium - 570 mm	<input type="checkbox"/> JAY030029 Tall - 610 mm

JAY J3 Back - Mid Deep Contour (150 mm Contour Depth)

Mid Thoracic:	<input type="checkbox"/> JAY030044	Medium - 340 mm	<input type="checkbox"/> JAY030041	Tall - 380 mm	
Upper Thoracic:	<input type="checkbox"/> JAY030045	Short - 420 mm	<input type="checkbox"/> JAY030042	Medium - 460 mm	<input type="checkbox"/> JAY030046 Tall - 500 mm
Shoulder Height:	<input type="checkbox"/> JAY030043	Short - 530 mm	<input type="checkbox"/> JAY030047	Medium - 570 mm	<input type="checkbox"/> JAY030048 Tall - 610 mm

Jay 3 Carbon Back - Shallow Contour

Lower Thoracic:	<input type="checkbox"/> JAY060001	Short - 170 mm	<input type="checkbox"/> JAY060002	Tall - 240 mm	
Mid Thoracic:	<input type="checkbox"/> JAY060003	Short - 300 mm	<input type="checkbox"/> JAY060004	Tall - 380 mm	

Accessories for JAY J3 backs

WHT030415	<input type="checkbox"/> Whitmyer Headrest, COBRA XTRA, flip back	(only with JAY J3 backs)	
WHT030416	<input type="checkbox"/> Whitmyer Headrest, AXYS, Standard	(only with JAY J3 backs)	
JAY090058	<input type="checkbox"/> JAY Accessory bag	(only with Jay 3 Carbon backs)	

ARMREST AND SIDEGUARDS

Armrest

XES040005	<input type="checkbox"/> Desk sideguard, flip-up, removable	
XES040001	<input type="checkbox"/> Desk sideguard, armpad (250 mm), fixed, flip-up, removable	
XES040002	<input type="checkbox"/> Desk sideguard, armpad (250 mm), height-adjustable, flip-up, removable	
XES040003	<input type="checkbox"/> Desk sideguard, armpad (310 mm), fixed, flip-up, removable	
XES040004	<input type="checkbox"/> Desk sideguard, armpad (310 mm), height-adjustable, flip-up, removable	
XES040052	<input type="checkbox"/> Single-post sideguard, armpad (330 mm), depth adjustable, height-adjustable	
XES040117	<input type="checkbox"/> Tubular armrest, Aluminium, swing away (sideways), padded, removable	

Accessories armrest

XES040020	<input type="checkbox"/> Release assistance for tetraplegics	(only for desk sideguard)	
XES090011	<input type="checkbox"/> Tray, push-fit, SW 380 - 480 mm		
	<input type="checkbox"/> Tray, swing-away (to the side), SW 340 - 500 mm		
XES090012	<input type="checkbox"/> left-mount	XES090013	<input type="checkbox"/> right-mount

Sideguards

XES040107	<input type="checkbox"/> Sideguard, composite, removable, 210 mm high (310 mm long), black	
XES040101	<input type="checkbox"/> Sideguard, carbon fibre, cold weather protection, slim (black)	
XES040102	<input type="checkbox"/> Sideguard, aluminium, cold weather protection (in frame colour)	
XES040103	<input type="checkbox"/> Sideguard, aluminium, fender, cold weather protection (in frame colour)	
XES040104	<input type="checkbox"/> Sideguard, carbon fibre, fender, cold weather protection (black)	

FOOTREST & FOOTBOARDS

XES050008	● Hanger, aluminium, swing-away, pair	(angle measured seat and hanger tube)	
XES050001	<input type="checkbox"/> Hanger angle 100°	XES050002	<input type="checkbox"/> Hanger angle 110°

Hanger, aluminium, flip-up (0° - 90°), swing-away

XES050009	<input type="checkbox"/> left	XES050010	<input type="checkbox"/> right
XES050011	<input type="checkbox"/> pair		

Amputee support

XES050014	<input type="checkbox"/> left	XES050013	<input type="checkbox"/> right
XES050012	<input type="checkbox"/> pair		

For more information
visit our webpage
www.SunriseMedical.eu
© Sunrise Medical 2023

JAY J3 Shallow Contour



JAY J3 Mid Contour




JAY J3 Mid Deep Contour



XES040107
Composite
sideguard
detachable



XES040102
Aluminium, cold weather
protection



XES040103
Aluminium, cold weather
protection, fender



XES040117
Tubular armrest, swing away,
removable, padded



XES04000x
Desk sideguard, armpad, flip-
up, removable



XES040104
Carbon fibre, cold weather
protection, fender

Xenon² SA

● = Standard ○ = No cost option □ = Cost option

Price in €

Footboards

XES050020	<input type="checkbox"/>	Footboard, divided, composit, angle-adjustable, flip-up (sideways), heel loop, black (from SW 380mm)	
XES050022	<input type="checkbox"/>	Footboard, divided, aluminium, angle-adjustable, flip-up (sideways), heel loop, frame colour	
XES050051	<input type="checkbox"/>	Footboard, divided, aluminium, angle-adj., flip-up (sideways), high-mount, heel loop, frame colour	
XES050023	<input type="checkbox"/>	Footboard, platform, aluminium, angle-adjustable, flip-up (sideways), heel loop, frame colour	
XES050126	<input type="checkbox"/>	Swing-away mechanism for platform footboard (minimum LLL 30mm longer)	
XES050059	<input type="checkbox"/>	Footboard, platform, Performance, aluminium, angle-adjustable, flip-up (sideways), calf strap	
XES050131	<input type="checkbox"/>	Footboard, platform, lightweight, plastic, angle & depth-adjustable flip up (sideways), calf strap, black	
XES050132	<input type="checkbox"/>	Footboard, platform, lightweight, carbon fibre, angle & depth-adj. flip up (sideways), calf strap, black	
XES050052	<input type="checkbox"/>	Footboard, platform, aluminium, angle-adjustable, flip-up (sideways), high-mount, frame colour	

Accessories for footboards

XES050127	<input type="checkbox"/>	Adjustable side protection, to create a side-stop on the footboard (1)	
-----------	--------------------------	--	--

LOWER LEG LENGTH (LLL): LLL has to be 25 mm less than front seat height; all measurements are in mm

<input type="radio"/>	220	<input type="radio"/>	230	<input type="radio"/>	240	<input type="radio"/>	250	<input type="radio"/>	260	<input type="radio"/>	270	<input type="radio"/>	280	<input type="radio"/>	290	<input type="radio"/>	300	<input type="radio"/>	310
<input type="radio"/>	320	<input type="radio"/>	330	<input type="radio"/>	340	<input type="radio"/>	350	<input type="radio"/>	360	<input type="radio"/>	370	<input type="radio"/>	380	<input type="radio"/>	390	<input type="radio"/>	400	<input type="radio"/>	410
<input type="radio"/>	420	<input type="radio"/>	430	<input type="radio"/>	440	<input type="radio"/>	450	<input type="radio"/>	460	<input type="radio"/>	470	<input type="radio"/>	480	<input type="radio"/>	490	<input type="radio"/>	500		

FRAME VERSIONS

XES010015	●	Frame SA (removable footrest)	
-----------	---	-------------------------------	--

CENTRE OF GRAVITY (COG)

XES070703	Adapter, standard	<input type="radio"/>	2	<input type="radio"/>	3	<input type="radio"/>	4	<input type="radio"/>	5	<input type="radio"/>	6		
		<input type="radio"/>	7	<input type="radio"/>	8	<input type="radio"/>	9	<input type="radio"/>	10				
XES070701	Adapter, lightweight (2)	<input type="radio"/>	5	<input type="radio"/>	6	<input type="radio"/>	7	<input type="radio"/>	9	<input type="radio"/>	10	<input type="radio"/>	11
XES010411	<input type="radio"/>	Axle plate, aluminium, long (full adjustment)											
XES010412	<input type="checkbox"/>	Axle plate, aluminium, shorted (max. adjustment 10 mm)											

FRAME ACCESSORIES

(choose either tip assist or anti- tip tubes per side)

XES090016	<input type="checkbox"/>	Frame protection, black	
XES090100	<input type="checkbox"/>	Tie down brackets for transportation	
XES090010	<input type="checkbox"/>	Transit wheels	
XES090050	<input type="checkbox"/>	Anti-tip tubes, plug-in, pair	
Anti-tip tubes, swing-away			
XES090004	<input type="checkbox"/>	Left XFF090005	<input type="checkbox"/> Right
XES090006	<input type="checkbox"/>	Pair	
Tip assist			
XES090001	<input type="checkbox"/>	Left XFF090002	<input type="checkbox"/> Right
XFF090003	<input type="checkbox"/>	Pair	
XES090009	<input type="checkbox"/>	Crutch holder, loop	(only in combination with tip assist)

EMPULSE F55 BUNDLING

EMO000101	<input type="radio"/>	With F55 order, not mounted to the wheelchair (Please use separate OF to order the F55.)	
EMO000102	<input type="radio"/>	With F55 order, mounted to the wheelchair (Please use separate OF to order the F55.)	

FORKS

(compare with your choice on the first page)

XES080029	<input type="radio"/>	Castor fork, Aluminium (only 138 & 174 mm fork)	
XES080013	<input type="radio"/>	Castor fork, Carbotecture® Colour: Black	
		Colour: <input type="checkbox"/> Silver <input type="checkbox"/> Orange <input type="checkbox"/> Red (3) <input type="checkbox"/> Blue (3)	
XES080008	<input type="checkbox"/>	Castor fork, aluminium, one-arm Colour: Black	
		Colour: <input type="checkbox"/> Silver <input type="checkbox"/> Orange <input type="checkbox"/> Red <input type="checkbox"/> Blue	
XES080007	<input type="checkbox"/>	Castor fork, aluminium, suspension (Frog Legs®), black	

- (1) only for footboard, platform lightweight or Performance
 (2) only with 0° wheel camber
 (3) only in combination with castor, soft (XES080302 and XES080306)

For more information
 visit our webpage
www.SunriseMedical.eu
 © Sunrise Medical 2023



XES050132
Footboard platform,
lightweight, carbon fibre



XES050059
Footboard platform,
Performance



XES050022
Footboard divided,
aluminium, frame colour



XES050127
Adjustable side protection



XES090001
Tip assist



XES090004
Anti-tip tubes, swing away



XES090050
Anti-tip tubes, plug in, pair



XES090010
Transit wheels

COG	Seat Depth						
	380	400	420	440	460	480	500
40	✓	✓	✓	✓	✓	✓	✓
60	✗	✓	✓	✓	✓	✓	✓
70	✗	✗	✓	✓	✓	✓	✓
80	✗	✗	✓	✓	✓	✓	✓
90	✗	✗	✗	✓	✓	✓	✓
100	✗	✗	✗	✓	✓	✓	✓
110	✗	✗	✗	✗	✓	✓	✓
120	✗	✗	✗	✗	✓	✓	✓
130	✗	✗	✗	✗	✗	✓	✓

Xenon² SA

● = Standard ○ = No cost option □ = Cost option

Price in €

CASTORS

(compare with your choice on the first page)

XES080305	○	Castor, pneumatic, rim: composite, black	(only 7")	
XES080301	○	Castor, solid (hard), rim: composite, black		
XES080302	○	Castor, solid (soft), rim: composite, black	(only 5" & 6")	
XES080303	□	Castor, solid (hard), light-up, rim: composite, black	(only 4" & 5")	
XES080304	□	Castor, solid (hard), rim: aluminium Colour, rim: Silver	(only 4" & 5")	
XES080306	□	Castor, solid (soft), rim: aluminium Colour, rim: Silver	(only 4" & 5")	
		Colour, rim: □ Black □ Orange □ Red □ Blue		

REAR WHEELS

(compare with your choice on the first page; choose only one wheel type)

XES070000	○	Without rear wheel		
XES070001	○	Rear wheel, Universal, 36 spokes (silver), cross-wise spokes	(only 24")	
XES070002	○	Rear wheel, Design, 36 spokes (silver), straight spokes	(only 24")	
XES070003	○	Rear wheel, Lightweight, 24 spokes (silver), straight spokes	(please see colour options on last page)	
XES070004	□	Rear wheel, Proton, 24 spokes (black), straight spokes	(user weight limit 100kg)	
XES070005	□	Rear wheel, Spinergy, 18 spokes (black), straight spokes; hub: silver		
XES070007	□	Rear wheel, Spinergy, 12 spokes (black), straight spokes, hub black		
XES070008	□	Rear wheel, drum brakes, 36 spokes (silver), cross-wise spokes (1)	(only 24")	
XES070009	□	Rear wheel, Mountainbike, 36 spokes (silver), straight spokes	(only 24")	
XES070010	□	Mountainbike wheel, replacement kit incl. tyre (extra coarse profile)	(only 24")	

Axle

XES070150	●	Quick-release axle		
XES070151	□	Quick-release axle, designed for quads (2)		

TYRES

(choose only one tyre type)

XES070101	○	Tyre, pneumatic, fine profile, Schwalbe Right Run		
XES070102	○	Tyre, puncture proof, tread		
XES070107	□	Tyre, pneumatic, slick, Schwalbe One		
XES070104	□	Tyre, pneumatic, puncture proof, tread, Schwalbe Marathon Plus Evolution	(only 24")	
XES090163	□	Spares-kit: 1 extra tyre (same as on the chair) and 1 extra inner tube	(not with solid tyres)	

HANDRIMS

XES070204	○	Handrims, stainless steel, silver		
XES070202	○	Handrims, aluminium, anodised, silver		
XES070201	○	Handrims, aluminium, powder-coated, clear		
XES070203	□	Handrims, aluminium, powder-coated, coloured	(select colour on last page)	
XES070205	□	Handrims, titanium, untreated		
XES070206	□	Handrims, titanium, powder-coated, clear		
XES070207	□	Handrims, Supergrip ®		
XES070208	□	Handrims, Maxgrepp ®	(only 24")	
XES070212	□	Handrims, Ellipse 3R, black	(ergonomic handrims with rubber-strip, 24" and 25" only)	

Hand rim mounting

XES070307	□	Handrim mount, super narrow (7 mm)	(only with stainless steel)	
XES070317	●	Handrim mount, narrow (17 mm)		
XES070327	○	Handrim mount, wide (27 mm)		

Handrim accessories

XES070250	□	Handrim cover, rubber	(only 24")	
-----------	---	-----------------------	------------	--

WHEEL CAMBER

(choose only one option)

XES070600	○	0°		
XES070602	○	2°	(not with Centre of Gravity Adapter, lightweight)	
XES070604	○	4°	(not with Centre of Gravity Adapter, lightweight)	

- (1) not with push handles, fold down & not without push handles.
only with Adapter, standard and only with 0° and 4° wheel camber
- (2) not with Rear wheels Spinergy/ Mountainbike/ drum brakes

For more information
visit our webpage
www.SunriseMedical.eu
© Sunrise Medical 21023



*XES070001
Universal wheel*



*XES070002
Design wheel*



*XES070003
Lightweight wheel*



*XES07000
Proton wheel*



*XES07010
Mountainbike wheel*



*XES070101
Schwalbe Right Run*



*XES070107
Schwalbe One*



*XES070104
Schwalbe Marathon Plus
Evolution*



*XES070102
Tyre, puncture proof*



*XES070202
Handrims aluminum,
anodized*



*XES070203
Handrims aluminum,
powder coated, coloured*



*XES070208
Handrims Maxgrepp*



*XES070212
Handrims Ellipse 3R*

Xenon² SA

● = Standard ○ = No cost option □ = Cost option

Price in €

SPOKE GUARDS

(only 24")

XES070831	<input type="checkbox"/>	Transparent		
		XFF070813	<input type="checkbox"/>	Galactica
XES070823	<input type="checkbox"/>	Princess	XFF070818	<input type="checkbox"/> Jungle Fever
XES070816	<input type="checkbox"/>	Happy Lions	XFF070811	<input type="checkbox"/> Feather
XES070819	<input type="checkbox"/>	Monkeys	XFF070821	<input type="checkbox"/> Pirates
XES070820	<input type="checkbox"/>	Oriental	XFF070806	<input type="checkbox"/> Chrome
XES070829	<input type="checkbox"/>	Superhero	XFF070814	<input type="checkbox"/> Graffiti
XES070817	<input type="checkbox"/>	Hot Wheels	XFF070802	<input type="checkbox"/> Basketball
			XFF070803	<input type="checkbox"/> Basketball 2
XES070822	<input type="checkbox"/>	Pool	XFF070812	<input type="checkbox"/> Flying Water
XES070827	<input type="checkbox"/>	Silk	XFF070828	<input type="checkbox"/> Space
XES070830	<input type="checkbox"/>	Tiger	XFF070833	<input type="checkbox"/> Zebra

WHEEL LOCKS

(choose only one option)

XES060001	●	Standard brake		
XES060002	<input type="checkbox"/>	Knee-lever brake		
XES060010	<input type="checkbox"/>	Brake lever extension, knee-lever brake, fold-down		
XES060023	<input type="checkbox"/>	Compact brake EVO		

ACCESSORIES

XES090000	●	Tool kit		
XES090025	<input type="checkbox"/>	Air pump		
XES090024	<input type="checkbox"/>	High pressure pump (0 - 11 bar)		
XES090030	<input type="checkbox"/>	Back pack		
XES090026	<input type="checkbox"/>	Mobile phone pocket		25 €

COLOURS ⁽¹⁾

Please consult our colour catalogue to review your colour choice

	Frame	Handrim	Hub/rim	X-brace		Frame	Handrim	Hub/rir	X-brace
02 Brilliant Silver	○	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	27 Graphitemetallic	○	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
03 Silver Matt	○	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	29 Black	○	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
06 Sun Yellow	○	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	31 Matt Black	○	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
06 Sun Yellow	○	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	32 Signal White	○	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
14 Blue	○	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	33 Pink	○	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
15 Brilliant Marine	○	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	34 Dark Chocolate Metallic	○	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
16 Matt Marine	○	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	35 Skyblue	○	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
19 Brilliant Purple	○	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	36 Cyan-Metallic	○	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
21 Red Ruby	○	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	37 Marsred (Orange)	○	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
22 Brilliant Red	○	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	40 Kiwi	○	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
23 Matt Red	○	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	41 Flamingo-Glitter	○	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
24 Tutti Frutti Orange	○	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	42 Matt-Anthracite	○	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
46 Mamba Green (new)	○	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>					
					Rear wheel, lightweight anodized, orange			<input type="checkbox"/>	
					Rear wheel, lightweight anodized, blue			<input type="checkbox"/>	
					Rear wheel, lightweight anodized, red			<input type="checkbox"/>	

COLOURS - for stickers ("Xenon²" on the frame and "Quickie" on sideguard)

XES100550	○	black (3)		-
XES100551	○	silver		-
XES100552	○	orange		-
XES100553	○	red		-
XES100554	○	blue		-

(1) Handrim: only with handrim colour powder coated; Hub/rim: only for universal wheel 24"

(2) only with 24" Lightweight wheel, Hub and rim anodized, black spokes

(3) Sticker "Quickie" on sideguard will be silver in combination with carbon sideguards

For more information
visit our webpage

www.SunriseMedical.eu

© Sunrise Medical 2023



*XES060001
Standard brake*



*XES060002
Knee lever brake*



*XES060023
Compact brake EVO*



*XES090000
Tool kit*



*XES090024
High pressure pump*



*XES090025
Air pump*



*XES090030
Back pack*



*XES090026
Mobile phone pocket*

Xenon² SA

● = Standard ○ = No cost option □ = Cost option

Price in €

Comments for special build - please apply for a quote; lead time conditions may change.



Quickie Xenon² SA - Basic Price

Option price	 €
Packaging and transport	 €
Exclusive VAT	 €
VAT.	 €
Total Price	 €

Sunrise Medical Customer service

Address: Kahlbachring 2-4, D-69254 Malsch/Heidelberg

Tel. +497253 980-0

Fax. +497253 980-222

E-mail: export@sunrisemedical.de

Website: sunrisemedical.eu

For more information

visit our webpage

www.SunriseMedical.eu

© Sunrise Medical 2023