



zippie®
Explore your world!



Youngster 3



Fit for Adventures



Youngster 3

The lightest folding wheelchair
with growth adjustability!

Zippie® Youngster 3™:

The new lightweight Youngster 3 offers integrated seat width and seat depth growth. The new frame design gives your wheelchair a stylish and unique look.

On top of that the Youngster 3 offers new designed, innovative and lightweight options using latest technology.





Patented 2 cm seat width growth via cross brace for growing user (without new parts).



New lightweight fork and frame for maximum adjustability. Improved look and better manoeuvrability.

- Lightweight package:

Lightweight fork, lightweight seat upholstery, lightweight anti-tips, Safari package



4 cm seat depth growth via seat tubes for growing user (without new parts).

Explore your world! Have fun, enjoy playing with your friends and be active

Options

The new Youngster 3™ offers you and your parents a wide range of options for more comfort, freedom, fun and safety.



New selection of back upholsteries with a wide range of different versions and an improved adjustability.



Lightweight seat upholstery for easy handling and transport.

Twentyeight frame colours!

Select a colour for each component. Cushions are always black




Modern sideguard design in combination with the new Safari wheel lock for more safety.

Spokeguards

There are 32 fantastic spokeguard designs you can choose from. Here you can see 4 designs! See our spoke guard catalogue for the whole range.































Technical specifications



Frame colour :	28 colours available (see order form)	Weight of heaviest part:	4.4 kg
Upholstery colour:	Black / Vented (black & silver) / Alcantara (grey & brown)	Product weight:	from 8.5 kg
Max user weight:	85 kg	Front wheel diameter:	4" / 5" / 6"
Growth: Seat width:	+2 cm without parts	Rear wheel diameter:	20" / 22" / 24"
Seat depth:	+4 cm without parts		

Frame Colours

							
1 White RAL 9010	2 Brilliant Silver	3 Silver matt	4 Silver	6 Sun Yellow	7 Lime	8 Matt Lime	9 Candy Teal
							
10 Matt Candy Teal	13 Flip Flop Blue/Green	14 Blue	15 Brilliant Marine	16 Matt Marine	19 Brilliant Purple	20 Matt Purple	21 Red Ruby
							
22 Brillant Red	23 Matt Red	24 Tutti Frutti Orange	26 Silver texture	27 Graphite-metallic	28 Platinmetallic	29 Black	30 Mystery
							
31 Matt Black	32 Signal white	33 Pink	34 Dark chocolate metallic				

It is best to select colours using the metal colour plates at your dealers and to use the order form to finalise your desired colour scheme



PB_Youngster3_EN_EU_REV.1

Sunrise Medical GmbH

Kahlbachring 2-4
69254 Malsch/HD
Tel.: +49 (0) 72 53 / 9 80-0
Fax: +49 (0) 72 53 / 9 80-111
www.SunriseMedical.eu