

GMFCS
Levels II–IV



Backrest
can be adjusted from 90°
to 180°.

Pelvic harness
provides reliable retention.

Seat surface
moulds to fit body
proportions.

I am your mobile bathing seat. I would really love to go bathing with you. But not only that: I would also be delighted to build castles with you in the sandbox. Hair full of sand, mucky hands – we have earned ourselves a shower! My backrest can be adjusted so that you can either sit bolt upright or lean back and relax. By the way: I am so light-weight and compact in size that you just have to fold me up. That's really practical when I accompany you on journeys or to friends' houses.



Back angle adjustment
brings your back into the
required position.



Data	bibi.	bibi. with long back
Seat depth	30 cm	30 cm
Seat width	36 cm	36 cm
High back	45 cm	55 cm
Total length	60 cm	60 cm
Back angle	90° - 180°	90° - 180°
max. load	30 kg	30 kg
Weight	2.3 kg	2.3 kg
HMV No.	04.40.04.1xxx	

Aluminium frame

It can be easily folded up to save space.



You can find more information on accessories & harnesses on the pages in the Price List.

