

On the same wavelength.

GMFCS
Levels III-V

kg
max. 114



I am your universal bathing and shower bench.

I can relax beautifully on my comfortable seating surface. Here you can adjust my back so that it is best for you – even directly whilst bathing. So that you also feel safe, there are different straps available for your chest, pelvis or legs. Of course, there is a headrest too. You can have me with or without a calf pad. My transfer sub-frame offers you ideal support when showering over the bathtub.

Bath sub-frame

puts you at the right height and relieves the strain on your parents' backs.



Shower sub-frame

with wide castors which can be braked; can also be used for transfer.





Headrest

can be adapted in height and width.

Removable cover

for fast and uncomplicated cleaning in the washing machine.

Plastic sub-frame

can be adjusted in height.

Support feet

can be adjusted in height.

**Chest strap**

provides security and support. This is available with and without the adjustable side guidance and can also be used as a pelvic strap.

Optional leg strap

is adapted appropriately using hook and eye tape.

Optional calf pad

can be adjusted during use.

Back angle

can be changed using one hand.



The angle of the calf pad and seat surface can be independently adjusted even during use – so that you can sit stably and comfortably or lie down.



Frame colour

○ White



Using the transfer frame underneath me, I can easily position you over the bathtub and even swivel by 90°. That makes showers such fun!



Data	Size 0	Size 1	Size 2
Seat depth	27 cm	33 cm	43 cm
Seat width	38 cm	38 cm	44 cm
Back height	58 cm	67 cm	81 cm
Seat height	5 cm / 18 cm	5 cm / 18 cm	5 cm / 18 cm
Seat height with bath sub-frame	32 - 47 cm	32 - 47 cm	32 - 47 cm
Seat height with shower sub-frame	70 cm	70 cm	70 cm
Seat height with transfer sub-frame	53 - 75 cm	53 - 75 cm	53 - 75 cm
Length of bathing bench, completely flat, without calf pad	89 cm	104 cm	128 cm
Length of bathing bench, completely flat, with calf pad	108 cm	124 cm	155 cm
Width	46 cm	46 cm	53 cm
Length, completely flat	75 cm	75 cm	75 cm
Height of bath sub-frame	20 cm	20 cm	20 cm
Length of bath sub-frame	99 cm	99 cm	99 cm
Width of bath sub-frame	44 cm	44 cm	44 cm
Height of shower sub-frame	56 cm	56 cm	56 cm
Width of shower sub-frame	69 cm	69 cm	69 cm
Length of shower sub-frame	86 cm	86 cm	86 cm
max. load	34 kg	68 kg	114 kg
Back and seat angle	0°, 15°, 30°, 45°, 60°, 75°		
HMV No.	04.40.04.1001		



You can find more information on accessories & harnesses
on the pages
in the Price List.

