

# siggi combi isofix.

## Buckle up, and off we go!

GMFCS  
Levels III-V

kg  
max. 18



20°–35°

25 cm

27 cm

new



### I am your special car seat, and stay reliably at your side during every ride in the car.

With me, you can really relax when underway, even if the journey takes a little longer. That's because my seat can be adjusted in four steps to fit you. In this way, you not only sit comfortably, but can also switch to a reclining position. The best thing about it is that you can

install me any way you like – either facing away from the

direction of travel using the Isofix (0–18 kg) or in the direction of travel without the Isofix (9–18 kg).

In this way, you and your parents can enjoy me for a long time to come. I like to boast, too – I have proven my skills in numerous independent tests!



0 to  
18 kg

semi-universal  
0–18 kg  
E5  
E5 04204

universal  
9–18 kg  
E5  
E5 04205

Special needs  
restraint 0–18 kg  
E5  
031 80S



# siggi combi isofix.

## Four-step

### angle adjustment

facilitates reliable and relaxed transport.

## 5-point strap

facilitate optimum hold.

## Isofix attachment points

for use facing away from the direction of travel.

## Headrest

can be adjusted in height.

## Thorax pelotte pads

can be adjusted in height and width.

## Front bar

to secure the child car seat and for more freedom to move your legs.

## Low tare weight

ensures easy installation and removal.

Child faces  
facing away  
from the direc-  
tion of travel.

0 to  
18 kg



Child faces  
in the direction  
of travel.

9 to  
18 kg





Facing away  
from the direc-  
tion of travel up  
to 18 kg.



Data		siggi combi isofix.
Seat width at front		32 cm
Seat width at rear		27 cm
Seat surface depth		25 cm
Height of headrest		48 - 60 cm
Seat angle	facing the direction of travel	30° / 37.5° / 45°
	facing away from the direction of travel*	25° / 32.5° / 40°
Total width		46 cm
Total height		68 cm
Group		0 - 18 kg
Weight		16 kg
HMV No.		26.11.06.0021

\* at a vehicle seat angle of 5°



Integrated  
magnets

hold the straps tight so that you  
can simply be placed into the seat.



You can find more infor-  
mation on accessories  
in the Price List.

