

# Simply outgrow yourself.

**GMFCS**  
Levels III-V



## I am the seat system which grows with you.

My dearest wish is to accompany you for a long time to come. For this reason, I am available in four different sizes which grow along with you. You can go anywhere with me, whatever the occasion - I am so flexible that I fit onto lots of different sub-frames. Apropos flexibility: you can adjust my pelvic guidance as you wish. Or simply remove my back unit.

### Headrest

can be adjusted in height, depth and angle.

### Back unit

optionally removable, can be adjusted in angle and height.

### Thorax pelotte pads

facilitate a stable and up-right sitting position.

### Armrests

can be adjusted in position.

### Thigh guide

supports the movement of the upper thighs outwards.

### Abducting leg guides

optionally available.



#### Frame colours



Lazurite blue  
metallic



Pure white



White aluminium



Anthracite metallic



Dormant red



With diverse  
accessories!



Data	Size 0	Size 1	Size 1b	Size 2
Seat depth	15 - 25 cm	21 - 30 cm	26 - 36 cm	30 - 42 cm
Seat width	16 - 26 cm	20 - 30 cm	20 - 30 cm	28 - 38 cm
Back height (short back)	35 - 43 cm	35 - 43 cm	42 - 52 cm	42 - 52 cm
Back height (long back)	42 - 52 cm	42 - 52 cm	53 - 63 cm	53 - 63 cm
Back angle	(-)5° - 25°	(-)5° - 25°	(-)5° - 25°	(-)5° - 25°
max. load	30 kg	30 kg	45 kg	50 kg
Weight	9 kg	10 kg	10 kg	12 kg
Recommended for sub-frame(s)	<b>freddy.</b> Size 0	<b>freddy.</b> Size 0	<b>freddy.</b> Size 1 + Anti-tip mechanism, <b>franka.</b> Size 2, <b>fridolin.</b> 1+2 Size 2	<b>freddy.</b> Size 1 + Anti-tip mechanism, <b>franka.</b> Size 2, <b>fridolin.</b> 1+2 Size 2
HMV No.	26.11.04.0xxx			



You can find more information on fabrics and embroidery  
on pages 174 - 175.

You can find more information on accessories and harnesses  
in the Price List.



# sofie & her friends.

## Growing together.



### fridolin 1.

the functional indoor sub-frame.



### fridolin 2.

The indoor sub-frame with height adjustment using a gas pressure spring.



### fridolin maxi.

the sturdy sub-frame for indoor use.



### freddy.

the compact seat pan sub-frame for outdoors.



### sofie.

the seat system that grows with the child.



### franka.

the sub-frame for indoor use with height adjustment using a gas pressure spring.

## As diverse as life itself.

Our seat system **sofie.** can be combined with very different sub-frames. In this way, you always have a seat tailored to your requirements for every situation or environment. A really clever system!



**sofie. &  
franka.**



**sofie. &  
fridolin 2.**



**sofie. &  
freddy.**