

**GMFCS**  
Level IV-V

  
kg  
max. 100

  
Home

  
(-)5° - 30°

  
44 - 56 cm

#### Support head

for all conventional seat pans and systems.


#### Foot pedal

facilitates inclination of the seat unit from (-)5° to 30°.

#### Dual castors

for easy transfer - with locking devices at the rear to prevent inadvertent rolling away.

**Central column**  
can be adjusted up to 12 cm in height.



**Height adjustment via adjustment wheel**



12 cm

**I am your robust sub-frame.** As my name implies, I am the biggest of us three brothers. That's why it is really easy for me to bring adolescents and adults weighing up to 100 kg from A to B. To get you up to the right height, my seat height can be adjusted using the adjustment wheel.

Frame colour



Anthracite metallic



## With a great range of accessories!

**Two gas pressure springs**  
of 200 N each can be operated  
by foot pedal



**Castors**  
with 100 mm  
width.



**Footrest**  
can be plugged on and is  
approx. 8 cm above the ground.



Data	fridolin maxi.
Height of seat pan support*	44-56 cm (46 - 58 cm**)
Seat angle	(-)5° - 30°
Length of base frame	80 cm
Total width	66 cm
Thickness of the gas pressure spring for tilt adjustment	2 x 200 N
max. load	100 kg
Weight	11,5 kg
HMV No.	26.99.01.0032

\* with 75 mm castors \*\* incl. trapezoidal adapter



You can find more information on accessories  
in the price list.

