



LIFE FAMILY

FIND THE ONE TO SUIT
YOUR LIFESTYLE

THE LIFE FAMILY

LIVE LIFE TO THE FULL.
THE EXTENDED LIFE FAMILY GIVES
YOU MORE FREEDOM THAN EVER.

The LIFE family impresses with its lightweight, robust, modern and highly adjustable wheelchairs.

With the choice of a rigid frame for performance or a folding frame for easier transfers, it is easy to find a LIFE to suit you.



LIFE

THE ALL-ROUNDER

Folding frame,
fixed & swing-away
front (shown)





LIFE R

THE SPORTY

Rigid frame,
fixed (shown) &
swing-away front



LIFE i

THE ADJUSTABLE

Folding frame,
swing-away front



A woman with short brown hair, wearing a light blue jacket and dark pants, is seated in a black wheelchair. She is positioned on a train platform, looking towards a blue and white train that has stopped with its doors open. The train's interior is visible through the windows, showing other passengers. The platform has a yellow safety line. The overall scene is brightly lit, suggesting an indoor or well-lit outdoor environment.

LIFE

THE WAY OF
MODERN LIFE

THE EASY LIFE

Live life to the full. Go out and explore the world.
Whatever your plans, your LIFE is ready.

Available as a fixed front or swing-away, the Quickie Life can cater for a number of needs and requirements in a foldable wheelchair. The Life utilizes a long established design that is proven to provide a very robust wheelchair and energy efficient set ups at a maximum user weight up to 140 kg.

But why stop there? Customise your LIFE by selecting from a wide range of side guards, armrests, hangers, footrests, wheels and handrims. Or define your style with a choice of 32 frame colours and anodized parts.



LIFE

THE ALL-ROUNDER

Folding frame,
Fixed & swing-away
front (shown)



Options and features of the LIFE – Configure this folding wheelchair to your individual needs. For further information on product options please check the LIFE order form.



// SEAT HEIGHT ADJUSTABLE IN SEVEN POSITIONS

With a camber from 0° to 4° and 4 centre of gravity positions (from the very active 8,3 cm to the more stable 2,3 cm).



// MODERN DESK ARMREST

Improved rigidity with a stylish design.



// FOLDING FIXATION

Operated with one hand, keeps the chair securely folded.



// MODERN SIDEGUARD DESIGN

The redesigned sideguards are available in aluminium or as stylish carbon version.



// ACTIVE SWING-AWAY FRONT

The swing-away hangers are now also available as 90° version for a more compact and sporty design.



// SMALL FOLDING PACKAGE

Options like the fold-down backrest or the removable swing-away hangers are leading to a very compact folding package.



// FLIP UP PLATFORM FOOTREST

Lightweight, angle and depth adjustable.



// INDIVIDUAL CUSTOMIZATION

Define your style with a choice of three different anodized style packages for the fork & axle plate.

LIFE R

THE PERFORMANCE FOR THE
SPECIAL MOMENTS IN LIFE.

THE ACTIVE LIFE R

Rigid, durable and robust, the LIFE R is designed to enhance your LIFE now and for years into the future.

Whatever the day throws at you, the energy efficient design of the LIFE R means you can focus on what's important to you.

Available as fixed front or swing-away, the new LIFE R has a modular design which means the front of the frame can be changed to increase seat depth, or alter the frame type.

The fixed front frame is available in different angles (75°, 85° and 85° with inset) and four frame lengths. The swing-away frame is available in three angles (70°, 80° and 90°) and a choice of hangers and three frame lengths.



LIFE R

THE SPORTY

Rigid frame,
Fixed (shown) &
swing-away front



Options and features of the LIFE R – Configure this rigid wheelchair to your individual needs. For further information on product options please check the LIFE R order form.



// SEAT HEIGHT ADJUSTABLE IN SEVEN POSITIONS

With a camber from 0° to 4° and four centre of gravity positions (from the very active 8,3 cm to the more steady 2,3 cm).



// UPDATED KNEE LEVER BRAKE

Easy to operate brake with swing away lever for easy transfers.



// ANGLE ADJUSTABLE BACK (OPTIONAL)

Lightweight and simple to adjust. Provides 20° of angle adjustment from -12° to +8° in 4° increments.



// ANODIZED 24" REAR WHEELS

Accentuate your colour scheme further by matching the rim and hub of your 24" rear wheel with 3 anodized colour selections.



// INDIVIDUAL CUSTOMIZATION

Define your style with a choice of three different anodized style packages for the fork & axle plate



// SURGE® LT HANDRIMS

Its Gripton rubber strip on top of the ergonomic oval shaped rim offers maximum traction.



// FOOTPLATE LOCK

This footboard remains always in position even when you are driving over rough terrain.



// ACTIVE SWING AWAY HANGERS

Offers a modern and active look. Available in 70°, 80° and 90° angle.



LIFE i

THE RIGHT ADJUSTMENTS
FOR INDIVIDUAL NEEDS.

THE ADAPTIVE LIFE i

Extensive adjustability makes LIFE i the ideal choice for users with progressive neurological disorders such as stroke, ALS or MS. Comprehensive adaptability ensures that the driving profile of the LIFE i is always in tune with changing user requirements.

Stroke patients are quite often affected by a hemiplegia. One arm drive and wheel lock option, half side tray and hemiplegic armrest are supportive options in this case. Also the low seat to floor height is intended for foot propelling users.

For MS (multiple sclerosis) and ALS (amyotrophic lateral sclerosis) patients our angle-adjustable and reclining backrest combined with the JAY backs options provide optimal trunk positioning. With its small increments of 1 cm for front and rear seat height the seat angle can be constantly adapted to changing needs.



LIFE i

THE ADJUSTABLE



Options and features of the LIFE i – Configure this rigid wheelchair to your individual needs. For further information on product options please check the LIFE R order form.



// EASY FRONT SEAT HEIGHT ADJUSTMENT

The full range of front seat height adjustment from 38 to 54 cm is designed for the neurological therapy. Ideal also for foot propelling users.



// MULTI-ADJUSTABLE REAR SEAT HEIGHT UNIT

With one unit the rear seat height, centre of gravity and camber can be easily adjusted.



// PHYSIOLOGICAL SEAT DEPTH GROWTH

The seat depth can be adjusted by moving the back tubes and the seat sling can be extended or shortened accordingly up to 4 cm.



// ONE-HAND OPERATION HEIGHT ADJUSTABLE ARMREST

The one-hand operation height adjustable armrest facilitates the height adjustment of the armpad especially for hemiplegic users.



// HALF SIDE TRAY

Half side tray e.g. for clients who suffered a stroke. Half Trays are available for the left or right side.



// HEMIPLEGIC ARMREST

This armrest provides good positioning of the arm especially for hemiplegic users.



// TRUNK SUPPORT

For higher trunk stability. The laterals are height, angle and depth adjustable.



// JAY BACKS AND CUSHIONS

The LIFE i can be equipped with the Jay products for advanced positioning.

TECHNICAL DATA

	Life	Life R	Life i
Max. user weight	140 kg	125 kg	140 kg
Frame version and front angle	Fixed front version: 75°/85 °/85° with 20 mm inset Swing-away front version: 90°/80°/70°/	Fixed front version: 75°/85 °/85° with 20 mm inset Swing-away front version: 90°/80°/70°/	Swing-away front version: 90°/80°/70°/
Seat width	360 - 540 mm	340 - 480 mm	340 - 540 mm
Seat depth	360 - 500 mm	360 - 500 mm	360 - 540 mm
Seat height front	400 - 540 mm	390 - 500 mm	380 - 540 mm
Seat height rear	400 - 510 mm	360 - 470 mm	370 - 510 mm
Backrest height	250 -500 mm	250 -450 mm	360 -500 mm
Backrest angle	78° - 102°	78° - 98°	85° - 115°
Wheel camber	0°/1°/2°/3°/4°	0°/1°/2°/3°/4°	0°/3°
Total width // length // height	760 mm // 1075 mm // 995 mm	720 mm // 1075 mm // 950 mm	800 mm // 1190 mm // 1125 mm
Max. pivot width	700 mm	700 mm	960 mm
Max. safe slope	7° (for use of parking brakes)	7° (for use of parking brakes)	7° (for use of parking brakes)
Front wheel // Rear wheel	4"/5"/6"/7" // 22"/24"/25"	4"/5"/6"/7" // 22"/24"/25"	5"/6"/7" // 22"/24"
Wheelchair weight	from 11 kg	from 11 kg	from 14,5 kg
Dismantled weight	from 6,7 kg	from 6,7 kg	from 10,2 kg
Mass of heaviest part	2,1 kg	2,1 kg	2,1 kg
Removable parts	Rear wheels & hangers (for transport only)	Rear wheels & hangers (for transport only)	Rear wheels & hangers (for transport only)
Colours	32	32	32
Crash tested	Yes (approved for transport in a vehicle according ISO 7176-19)		
Transport without user	Fix the wheelchair safe in the vehicle (e.g. car or plane)		
Intended use and environment	For people with limited mobility (children & adults); self-propelling or being pushed by an assistant, indoor and outdoor use.		



Sunrise Medical GmbH

Kahlbachring 2-4

69254 Malsch/HD

Tel.: +49 (0) 72 53 / 9 80-0

Fax: +49 (0) 72 53 / 9 80-222

E-Mail: international@sunrisemedical.de

www.SunriseMedical.eu

For further information on the full specification, options and accessories, please refer to the order form. All information is subject to change without notification. Please consult Sunrise Medical with any queries you may have. If you are visually impaired, this document can be viewed in PDF format at www.SunriseMedical.eu.

