



JAY®

JAY®



JAY EASY BACK

Simple, Comfortable and Effective

Choose your level of support

The JAY EASY Back system supports an upright posture and provides a stable and ergonomic seating position. The JAY EASY Back offers two different types of contours for an increased trunk stability depending on users needs.



Different conditions and needs of the client demand adequate solutions. The three different heights focus on users with a need for moderate to full posterior support.



JAY EASY BACK	
Backrest width	36 / 41 / 46 / 51 cm
Backrest height	38 / 46 / 53 cm
Back angle	+12° to -12°
Transport weight	2,4 kg
Max. user weight	136 kg
Accessories	Headrest, fixed lateral support, chest straps and lumbar supports
Safety	Crash tested



Sunrise Medical GmbH

Kahlbachring 2-4
69254 Malsch/HD
Tel.: +49 (0) 72 53 / 9 80-400
Fax: +49 (0) 72 53 / 9 80-222
E-Mail: international@sunrisemedical.de
www.SunriseMedical.eu



EASY to prescribe, EASY to order, EASY to use

With the JAY experience in seating and positioning, the JAY EASY backrest is designed to offer simple, comfortable and effective support.

EASY to prescribe: Just select from the two shapes and the three heights the right level of support. With the lumbar pads you can customize it further.

EASY to order: The hardware is designed to fit all common chair tubes. Width, depth and height adjustments are quick to do with just one tool.

EASY to use: The shell design combined with the microclimatic cover and the high quality foam provides the user with an effective and comfortable support that increases the seating tolerance.



Versatile One-step Release Hardware

The hardware is easy to operate and to adjust for an optimal back support.

- Backrest can be released with only one hand. Unlock system on both sides with the cord then angle backrest and pull up to take it out of the fixation
- Hardware fits on all commonly used back tubes: 1,9 cm with insert, 2,2 – 2,3 cm with insert or 2,5 cm without insert
- Adjustable in height, width, depth and back angle with just one tool



Innovative Cover Technology

JAY EASY Back comes with a microclimatic cover and a soft foam underneath for maximum comfort and user protection.

- 3DX™ spacer fabric vents both heat and moisture and increases air flow
- Cover is machine washable at 60° C and quick to dry
- Shape of the back can be easily adapted to the user through lumbar pads

++



Make it individual

Accessories – The JAY EASY Back comes with a useful range of accessories to give your back an optimal support.

- Headrest
- Lateral support
- Chest straps and lumbar supports

