

Weight capacity: 300 kg / 661 lbs

## Dimensions

Height 2300mm/91"

Length 2650mm/104" (71806)

Length 3150mm/124" (71808)

Width 998mm/39"

## Material

Rail, powder coated extruded aluminium (White)

Pole and foot, powder coated steel (White)

## Included parts:

- Rail 2,5 m or 3 m - 1pc
- Pole package, 2 pc
- Legs package, 2 pc
- Plate with screw 4pc
- Locking nut M10 4pc
- Locking plate 6pc
- Adjustable damper 2pc
- Rail roller 1pc
- Distance strap 1pc
- Rail nut 2pc
- Screw M10x30 2pc

## Tools for assembly

- Allen key 8mm and 6mm
- Philips screwdriver 5.5mm
- Flexible ended wrench 10 mm

## Product label

The product label contains information concerning the equipment.

The freestanding frame is tested by an accredited testing institute and comply with the requirements of directives for medical technical Class I products.



## Symbols on the unit



Type B product



Warning



User manual: Read and understand the manual before the product is used.



**ORTHOSTATICAL**  
*Live with autonomy!*