

TECHNICAL DATA

	LIFE T	LIFE RT
Max. user weight:	100 kg	100 kg
Seat width:	28 – 50 cm	28 – 48 cm
Seat depth:	32 – 46 cm	32 – 46 cm
Seat height front:	40 – 53 cm	36 – 50 cm
Seat height rear:	40 – 50 cm	36 – 45 cm
Backrest height:	25 – 47,5 cm	25 – 45 cm
Backrest angle:	78° to 102° in 4° increments	78° to 98° in 4° increments
Centre of gravity (COG):	2,3 – 8,3 cm	2,3 – 8,3 cm
Camber:	0° – 6°	0° – 6°
Fixed frame front:	75°; 75° with abduction; 85°; 85° with inset	75°; 75° with abduction; 85°; 85° with inset
Swing away frame front:	80° and 70°	80° and 70°
Chair weight:	from 10 kg	from 10 kg
Frame colours:	32	32
Crash tested:	yes	yes



For further information on product options please check the order form on our website.
You can view this document also in a larger font on a PDF-file at www.SunriseMedical.eu



Sunrise Medical GmbH

Kahlbachring 2-4
D-69254 Malsch/HD
Tel.: +49 (0) 72 53 / 9 80-0
Fax: +49 (0) 72 53 / 9 80-222
E-Mail: international@sunrisemedical.eu
www.SunriseMedical.eu

The background of the entire advertisement is a photograph of three teenagers. In the foreground, a young woman with blonde hair in a ponytail, wearing a light blue denim shirt and jeans, is seated in a blue QUICKIE wheelchair. She is smiling and looking towards the right. Behind her, another young woman with long brown hair, wearing a white t-shirt, is sitting on a red QUICKIE scooter and holding a tablet. In the background, a young man in a black t-shirt and jeans is standing and holding two wooden sticks, possibly for a game. The setting is outdoors, with a light blue corrugated metal wall on the left and green foliage in the background.

LIFE T & LIFE RT

MATCHING A TEENAGER'S ACTIVE LIFE



COOL, STYLISH, MODERN.

TWO WHEELCHAIRS – ONE FAMILY

Live life to the full! Meeting your friends at the mall? Having a fun time at the park? Whatever your plans are for the day, the LIFE RT and LIFE T are two wheelchairs that perfectly match the needs of a teenager's active life.

The LIFE RT is strong, reliable and offers the support you need to enjoy all kinds of cool activities. Need a folding frame? The LIFE T is not only easy to transport, but is also the perfect all-rounder for your daily routine. Regardless of your selection, both wheelchairs are designed for an adolescent's body size and needs. So whether you require a fixed-front frame (with a choice of straight, abduction or inset front ends) or you prefer a swing-away hanger version – you're sure to find a LIFE to suit you.



LIFE T
THE FOLDABLE

Folding frame,
fixed & swing away front



LIFE RT
THE SPORTY

Rigid frame,
fixed & swing away front



→ **Abduction frame for Life T & Life RT**

This frame version offers 2,5 cm abduction on each side. The abduction already starts below the seat upholstery.



LIFE T & LIFE RT – AN EXPRESSION OF A TEENAGER'S INDIVIDUALITY

Customise your wheelchair to your personal taste and needs. Define your style with a choice of more than 30 frame colours, three different style packages for the fork & axle plate and the full range of stylish spoke guards. Select from a wide range of side guards, armrests, hangers and footrests. You can even fine-tune your wheelchair with a sporty camber, stylish wheels and handrims or stunning Carbon options.

→ **Wide camber range for Life T & Life RT**

Camber improves the lateral stability of the wheelchair, protects the fingers when passing through narrow doorways and gives a sporty look. Both wheelchairs offer a camber range of 6° in 1° steps.





→ **The innovative Carbotecture fork**

Made out of stiff hightech material and extremely lightweight.
Available as a style package in three different colours.



→ **Carbon platform footplate**

choose a carbon footplate and a carbon sideguard to express your individuality.

LEONARDO'S KID



→ LIFE RT – the sporty



→ LIFE T – the foldable

Options and features of the LIFE T and LIFE RT – Configure your wheelchair to your individual needs. For further information on product options please check the LIFE T and LIFE RT order form.

→ For LIFE T & LIFE RT



// ABDUCTION FRAME

This frame version offers 2,5 cm abduction on each side. The abduction already starts below the seat upholstery.



// CARBON PLATFORM FOOTPLATE

Lightweight and stylish footplate with angle and depth adjustment.



// THE CONVINCING SAFARI WHEEL LOCK

With its intuitive handling, the safari sideguard wheel lock helps you stop quickly and safely.



// SEAT HEIGHT ADJUSTABLE IN SEVEN POSITIONS

With a camber from 0° to 6° and 4 centre of gravity positions (from the very active 8,3 cm to the more stable 2,3 cm).

→ For LIFE RT



// ANGLE ADJUSTABLE BACK

This angle adjustable backrest provides an adjustment from 78° to 98° in 4° increments.



// IT'S ALL ABOUT ERGONOMICS

The push bar can be easily adjusted in angle and height for an enhanced ergonomic position.



// DEPTH ADJUSTABLE SEAT UPHOLSTERY

Depth adjustable seat upholstery offers growth in combination with a long frame.

→ For LIFE T



// ANGLE ADJUSTABLE BACK

Lightweight and simple to adjust. Provides angle adjustment from 78° to 102° in 4° increments.