



75 - 95 cm



75 cm



58 cm



50 cm



34 - 28 / 28 - 22 cm



36 - 42 cm



6°



47 - 53 cm



66 cm



22 - 34 cm



90° - 112°



4 cm



75 kg



10 - 11 kg



Seduta da 28 riducibile fino a 22 per utenti da 5 a 3 anni circa.  
Seduta da 34 riducibile fino a 28 per utenti da 10 a 6 anni circa.

Seating 28 cm reducible until 22 cm for users from 6 to 3 years.  
Seating 34 cm reducible until 28 cm for users from 10 to 6 years.

## 17.56N BABY EVOLUTION

"Baby evolution" è una carrozzina molto leggera con carichi dinamici attivi. Con il ns. sistema 3D di variazione delle dimensioni da cm 28 adattabile da 3 a 5 anni circa e cm 34 adattabile da 10 a 6 anni circa, quindi è un prodotto che si adatta alle necessità di crescita degli utilizzatori. I fondi del sedile e dello schienale sono tensionabili con cuscino imbottito di appoggio; inoltre lo schienale può variare l'angolo fino a 12° ed è pieghevole in avanti. Le ruote grandi Ø 560 sono ad estrazione rapida; il telaio è in duralluminio; le ruote anteriori sono impermeforabili Ø 120 mm, le maniglie per l'accompagnatore sono regolabili ed asportabili; freni di stazionamento.

"Baby evolution" is a very light wheelchair with active dynamic loads. By our system 3D of adjustment of the dimensions from cm 28 adjustable from 3 to 5 years old and 34 cm adjustable from 10 to 6 years old, therefore it is a product which adapts to the growing needs of the users. The bottoms of the seating and of the backrest are composed by tensioning belts with laying padding cushion in addition the backrest can change the angle until 12° and is foldable forwards. The rear wheels Ø 560 are of the quick release type; the frame is made of duraluminium, the front castors are puncture proof type Ø 120, the handles for the attendant are adjustable and removable, parking brakes.



cm 28/22

cm 34/28

Particolari per regolazione e riduzione della seduta fino a -4 cm.

Detail for adjustment and seating reduction until -4 cm.



Particolari per regolazione fianchetti fino a -4cm

Detail for adjustment side panels until -4 cm.