

GMFCS
Levels III–V



I am your versatile bathing and shower bench. I will adapt completely to your requirements, whether in the bathroom, in the sea or on the beach. You will always lie safely and comfortably on me! My cover is washable to 30 °C, is particularly skin-friendly and leaves no pressure marks on your body. Because you can sit very low in the bathtub with me, I can save a lot of water. Aren't I practical? What's also practical is that my back angle can be adjusted with one hand during use – leaving the other hand free to support you!



Thanks to the head section, the angle of which can be changed, no water will run into your eyes.



Head section

can be easily folded.

Head guide pad

can be adjusted in height and width.

Cover

can be removed for cleaning purposes.

Backrest

can be adjusted from 90° to 180°.

Chest strap

provides reliable retention.

Knee angle adjustment

provides a comfortable sitting position.

Optional pelvic mould depression

stabilises and positioned the pelvis, with adjustable seat support for ischiatic bone.

Lower leg support

can be adjusted between 0° and 90°.

Sub-frame


for the shower.

**Pelvic mould depression**

optionally, a cover is available with an integrated pelvic mould depression.

Back angle

can be adjusted even during use with one hand.



In case of the Size 4, this takes place via two stainless steel snap joints.





Data	Size 1	Size 2	Size 3	Size 4
Seat width	30 cm	35 cm	35 cm	40 cm
Seat depth	30 cm	34 cm	38 cm	44 cm
Back height	45 cm	61 cm	71 cm	82 cm
Standing space length	69 cm	74 cm	88 cm	93 cm
Lower leg length	19 cm	20 cm	28 cm	28 cm
Back angle	90° - 180°	90° - 180°	90° - 180°	90° - 180°
Knee angle	0° - 90°	0° - 90°	0° - 90°	0° - 90°
Headrest angle	0° - 90°	0° - 90°	0° - 90°	0° - 90°
Headrest angle	0° / 5 cm	0° / 5 cm	0° / 5 cm	0° / 5 cm
max. load	35 kg	45 kg	55 kg	75 kg
Weight	4 kg	4,1 kg	4,3 kg	4,8 kg
HMV No.	04.40.04.1008			

* In combination with the cover with pelvic mould depression



You can find more information on accessories & harnesses
on the pages
in the Price List.

