

GMFCS
Levels III-V



I am your pivotable car seat. I love taking a ride in cars. Do you too? Then please take me with you! With me, you are certain to arrive safely. I have been tested and found to be absolutely safe - on the passenger seat and on the rear seat. To make it very easy for your parents to position you on my seat, it is child's play to pivot me by 90 degrees out of the car. And my best feature of all is that I can be used as the basis for special constructions. In this way, I can be equipped completely according to your requirements. There aren't many other systems around which can do that.

With a great
range of
accessories!





Seat cover

can be removed.

Thorax pelotte pads

facilitate a stable and upright sitting position.

Back

can be tilted in three steps.

Swivel device

can be pivoted by 90°.

Padded headrest

can be adjusted in height.

Straps

help to support and position.

Anatomic side pads

keep the body upright.

Foot support

is adjustable in height and can be removed for transfer.

90°

swivel device facilitates transfer.

Removable pelotte pads

Frame colour



Black – satin finished



Data	Size 1		Size 2		Size 3	
Seat depth	25 cm	30 cm	30 cm	35 cm	35 cm	40 cm
Seat width	28 cm		32 cm		36 cm	
Back height	40 cm		48 cm		56 cm	
Headrest height	59 - 69 cm		65 - 77 cm		76 - 86 cm	
Pelotte pad height	25 - 33 cm		32 - 40 cm		36 - 44 cm	
Seat angle*	0°/15°		0°/10°		0°/10°	
Chest width	14 - 25 cm		16 - 29 cm		16 - 32 cm	
Back angle	90° - 108°		90° - 108°		90° - 108°	
max. load	30 kg		45 kg		60 kg	
HMV No.	26.11.06.0009					

* When using the foot support and tilt adjustment



You can find more information on accessories in the Price List.

