

GMFCS  
Levels III–V

kg  
max. 120



(-)5°–25°

(-)5°–30°

36–47 cm

48–60 cm

34–44 cm



### **I am your sturdy therapy chair.**

I am your daily companion for adolescents and adults. Whether during your midday meal; when doing schoolwork or in reading hour: I will make sure you are in an optimum sitting position. Because my central column and my extended wheelbase guarantee a high level of stability, I can cope well with high loads. My anatomically-moulded seat surface ensures an upright pelvic position and makes sure that you can't slip off the seat surface.





# madita maxi.

**Push handle**  
for easy transfer.

**Thorax pelotte pads**  
facilitate a stable and up-  
right sitting position.

**Armrests**  
are equipped with extra-wide sup-  
ports and can be adjusted in width,  
height and angle.

**Foot pedal**  
for triggering the  
seat angle adjustment.

**Dual castors**  
facilitate transfer.

**Headrest**  
can be adjusted in height.

**Back angle**  
can be adjusted using a clam-  
ping lever.

**Thigh guide**  
can be adjusted in width.

**Abduction wedge**  
guides the legs.

**Sub-frame**  
is sturdy and robust.

**Simple, practical, good:**  
You can conveniently place  
your feet on the plug-in foot  
support, or independently jog  
along with the stool  
when the foot support is  
not installed.

**Easy  
assembly!**

Frame colour



Anthracite metallic



Data	Size 0
Seat depth	36 - 47 cm
Seat width	34 - 44 cm
Height of short back (made of 3 mm / 4 mm aluminium plate)	40 cm / 48 cm
Height of long back (made of 3 mm / 4 mm aluminium plate)	48 cm / 56 cm
Seat height	48 - 60 cm
Lower leg length	18 - 45 cm
Total width	64 cm
Length of base frame	80 cm
Seat angle	(-) 5° - 30°
Back angle	(-) 5° - 25°
max. load	120 kg
Weight	26 kg
HMV No.	26.11.05.0021

**Correct sitting made easy: your pelvis is supported thanks to the anatomically-moulded seat and back surface.**



**Anatomic seat and back surface.**



**You can find more information on fabrics and embroidery on pages 174 - 175.**

**You can find more information on accessories and harnesses in the Price List.**

