

madita-fun mini.

Small, but impressive.

GMFCS
Levels III-V


max. 20




(-)5° - 25°


(-)5° - 42°


12 - 19 cm


25 - 62 cm


14 - 22 cm

I am the therapy chair for when you are young, and give you security from the start.

I will place you on a level with your friends. My comfortable seat can be easily adjusted in height so that you can sit at any table and then be lowered back down to a height suitable for play. My optionally-available dysplasia seat permits you to sit correctly, so that your hip joints can develop positively. In addition, I offer you a narrow or wide back shape.





madita-fun mini.

Push handle

facilitates transfer.

Armrests

can be adapted in height and width.

Seat depth

extends from 12 to 19 cm.

Headrest

is available in a three-dimensionally adjustable or height adjustable version.

Thorax pelotte pads

facilitate a stable and up-right sitting position.

Seat

ensures optimum development of the hip joint.

Foot support

padded.

Seat height

can be adapted with support from a gas pressure spring.


Castors

facilitate transfer.

As individual as you are: in addition to a standard seat pad, you can select a pad with anatomic leg guidance or a dysplasia seat.



**Ideal for initial
care.**



Data	Size 0
Seat depth	12 - 19 cm
Seat width	14 - 22 cm
Height, short back	28 - 33 cm
Height, long back	30 - 35 cm
Seat height	25 - 62 cm
Lower leg length	12 - 22 cm
Total width	52 cm
Length of base frame	58 cm
Seat angle	(-)5° - 37° / 0° - 42°
Back angle	(-)5° - 25°
max. load	20 kg
Weight	16 kg
HMV No.	26.11.05.0022



**You can find more information on fabrics and embroidery
on pages 174 - 175.**

**You can find more information on accessories and harnesses
in the Price List.**

