

GMFCS
Levels III-IV

kg
max. 100



46-124 cm

29-113 cm



I am your versatile posterior walking

aid. I am always there for you – even from your first year of life. Walking is my biggest hobby. It doesn't matter where you want to go; I will follow you.

If you support yourself on my handles, you will quickly notice how I take on your body weight. In addition, my special abducted frame shape supports you so that you don't swing outwards when walking. My pelvic pelotte pad ensures that you walk upright and straight. Cool, aren't I?!





Buttocks pelotte pad
can become a backrest in
one single easy step.

Aluminium frame
can easily be folded
up.

Backwards stop
prevents inadvertent
rolling back.

EVA guide wheels
make the walking aid manoeuvrable.

Forearm rests
with handles.

Padded bench
can easily be folded
back.

Direction locks
facilitate
walking straight.



Abducted frame shape
supports you when walking
and prevents you from swinging
outwards.



Aluminium frame
can easily be folded up via
folding joints.

Frame colours
Size 0 mini



Frame colours
Size 1



Frame colours
Size 2



Frame colours
Size 3 and Size 4





Data	Size 0 mini*	Size 1	Size 2	Size 3	Size 4
Handle height	29 – 53 cm	34 – 58 cm	51 – 75 cm	68 – 94 cm	87 – 113 cm
Height of forearm rests	-	46 – 70 cm	63 – 87 cm	79 – 103 cm	98 – 124 cm
Handle spacing	22 – 36 cm	22 – 36 cm	30 – 44 cm	38 – 52 cm	43 – 57 cm
Seat height	-	33 cm	45 cm	60 cm	70 cm
Total width	52 cm	51 cm	61 cm	68 cm	75 cm
Total length	58 cm	58 cm	70 cm	77 cm	87 cm
Frame height (folded)	37 cm	25 cm	32 cm	34 cm	37 cm
Transport dimensions (L x W x H)	58 x 52 x 37 cm	58 x 51 x 25 cm	70 x 61 x 32 cm	77 x 68 x 34 cm	87 x 75 x 37 cm
max. load	20 kg	30 kg	50 kg	75 kg / 70 kg	100 kg
Weight	3.9 kg	7.1 kg	8.3 kg	8.5 kg / 9.8 kg	10.2 kg
HMV No.	10.46.02.3029				

* Frame cannot be folded up



You can find more information on accessories in the Price List.

