

momo motion.

For extra support.

GMFCS
Levels I-V

kg
max. 120



42-98.5 cm

18-23 cm

I am your tried-and-tested tricycle with E-drive.

I am here to help when you start to get tired. Whether on short or long trips into town or the countryside. Because I have a silent direct drive, you can enjoy your days out without being disturbed. In addition, my optional start-up assistance supports you, so that you get the tour started. For experienced drivers, my drive is also available in an "advanced version". This gives you a more powerful drive. Even my "beginner version" can be controlled in the Eco, standard or power mode using the display.





momo motion.



DirectPower E-bike motor

is designed brushless and without gears, minimises maintenance work and prevents motor noise.



LCD display

shows you information on the support level, speed and battery status - and is easy to read, even in bright / sunny surroundings.



Button or twist handle operation

allow you to set the mode, and the start-up assistance if necessary, including during travel.



Frame colours 16" & 20"

- Melon yellow
- Navy blue
- Pure white

Frame colours 24", 26" & 26" XL

- Anthracite metallic
- Pure white



**Wherever your travels take you, you can
always rely on me.
We are strong together!**



**The battery can also be
removed for charging.
The charge level can be read
on the display.**



	16"	20"	24"	26"	26" (XL)
Pedal to top edge of saddle*	42 – 57.5 cm	50 – 64.5 cm	56.5 – 73.5 cm	66.5 – 83 cm	74 – 90.5 cm
Pedal to top edge of saddle**	45 – 60 cm	52 – 68 cm	58 – 81.5 cm	67.5 – 91 cm	75 – 98.5 cm
Saddle support tube to front tube	39 cm	45 cm	52 cm	58 cm	58 cm
Wheelbase	80 cm	91 cm	104 cm	112 cm	112 cm
Access point height	18 cm	20 cm	22 cm	23 cm	23 cm
Full length	120 cm	140 cm	165 cm	178 cm	178 cm
Full width	65 cm	69 cm	75 cm	75 cm	75 cm
Max. load	60 kg	80 kg	100 kg	120 kg	120 kg
Weight	24.5 kg	27 kg	29.5 kg	31.5 kg	31.7 kg
HMV No.	22.51.04.0XXX (under application)				

* Dimensions with standard saddle post / ** Dimensions with T-saddle post



**You can find more information on accessories and harnesses
in the Catalogue and Price List.**

