

momo therapy bicycle.

Firmly in the saddle.

GMFCS
Levels I-V

kg
max. 100



30-95 cm

16-23 cm

I am your therapy bicycle for more

balance. I am there for you if the pleasant weather tempts you to take a trip around the lake. But I can do much more than that if you take me on an excursion. Depending on how steady you are on your bike, my stabilizer wheels can be adjusted from rigid to very soft spring settings. In this way, you can easily train your balance. Apart from that, I look really cool: without my supporting functions, I look just like a bicycle!





momo therapy bicycle.

Handlebars

Available in many different designs.

Individual drive

from a rigid drive via a 7-gear freewheel brake hub, right up to electrical support.

Saddle

available in different shapes and sizes.

Two brake systems

from 16", consisting of rims and drum brakes.

Pedals

Can be individually configured.

Stabilizer wheels

can be adjusted from rigid to very soft.

Low access level as standard feature

permits the rider to climb on conveniently and also allows independent transfer.

You can find a diverse range of accessories

for optimum support & individual equipment from page 122.

You have the choice! Between many handlebars, saddles, foot pans, head supports and different drives - from rigid toothed wheels right up to the 3-gear hub gear with reverse option!



Frame colours 12"

- Telemagenta
- Sky blue
- Pure white

Frame colours 16" and 20"

- Melon yellow
- Navy blue
- Pure white

Frame colours 24" and 26"

- Anthracite metallic
- Pure white



Using me, you can train your balance, coordination and endurance.



Dynamic stabilizer wheels

can be adjusted to your capabilities. If set to rigid, they will initially give you security and build trust. Over the course of your therapy, they can be adjusted to become very gentle suspension stabilizers. These can easily balance you as you go round bends and also offer sufficient drive force on dirt roads.



| Data | 12" | 16" | 20" | 24" | 26" |
|-----------------------------------|---------------|--------------|--------------|--------------|----------------|
| Pedal to top edge of saddle* | 30 - 38 cm | 42 - 57.5 cm | 50 - 64.5 cm | 59 - 76.5 cm | 71.5 - 87.5 cm |
| Pedal to top edge of saddle** | 33 - 41 cm | 45 - 60 cm | 52 - 68 cm | 61 - 84 cm | 73 - 95 cm |
| Saddle support tube to front tube | 33 cm | 39 cm | 45 cm | 53 cm | 58 cm |
| Wheelbase | 64 cm | 76 cm | 89 cm | 104 cm | 112 cm |
| Full length | 99 cm | 122 cm | 148 cm | 168 cm | 178 cm |
| Access point height | 16 cm | 18 cm | 20 cm | 22 cm | 23 cm |
| Full width | 43 cm | 52 cm | 62 cm | 71 cm | 76 cm |
| max. load | 35 kg | 50 kg | 65 kg | 100 kg | 100 kg |
| Weight | 10.5 kg | 13 kg | 15 kg | 18 kg | 20 kg |
| HMV No. | 22.51.01.0014 | | | | |



You can find further information on accessories & harnesses
on pages 122 – 125 and in the Price List.

