

freddy.

Wild and manoeuvrable.

GMFCS
Levels IV-V

kg
max. 120



23–50 cm

(-)5°–30°

25–35 cm

I am your compact seat pan sub-frame. You can go anywhere with me, whether to school or to friends. I am really manoeuvrable and very lightweight. You can adjust my height precisely according to your wishes. If you like, you can rest your feet on my adjustable supports. The consciously-designed short wheelbase creates maximum space for footrests and for your legs. And do you know what else is good about me? I have been tested in accordance with ISO 7176: Part 19 and therefore can also accompany you on car rides.



Individual footrests
with calf plates

can be adjusted in height, width and knee angle, and can also be swivelled away



Divided footrests

are adjustable in height and width and can be swivelled away.





**The scale applied to
the inner tubes**
allows you to determine
the set seating height.

Back
can be folded down onto the seat frame.

Spokes and clothing guard
prevent hands or clothing getting into the spokes.

Rear wheels
available in 16", 20", 22" and 24".

Push handle
for easy transfer.

Back angle adjustment
takes place via a snap joint or through central quick adjustment.

Scale on the inner tubes
for height adjustment of the seat adaptor.

Trapezoidal adapter
for all conventional seat pans.

Divided footrests
can be adjusted in height and in angle.

Tilting device
makes it easy to overcome kerbs.

Swing-away anti-tip mechanism
prevents you from tipping backwards.

Frame colours

- Lazurite blue metallic
- Pure white
- White aluminium
- Anthracite metallic
- Dormant red
- Black - satin finished



Data	Size 0	Size 1	Size 2	Size 3
min. seat height*	43 cm	43 cm	43 cm	43 cm
Short seat depth	23 - 33 cm	27 - 37 cm	-	-
Standard seat depth	28 - 38 cm	32 - 42 cm	40 - 50 cm	40 - 50 cm
Seat angle	(-)5° - 30°	(-)5° - 30°	(-)5° - 30°	(-)5° - 30°
Width between rear wheels (with clothing guard) = max. outer width of seat pan	39 cm	44 cm	49 cm	54 cm
Frame width inside	28 cm	33 cm	38 cm	43 cm
Frame width outside	35 cm	40 cm	45 cm	50 cm
Seat frame width	25 cm	30 cm	35 cm	35 cm
Total length**	63 - 79 cm	63 - 79 cm	66 - 82 cm	71 - 87 cm
Total width	58 cm	63 cm	68 cm	73 cm
max. load	90 kg	90 kg	120 kg	120 kg
Weight	From 14 kg	From 14 kg	From 14 kg	From 14 kg
HMV No.	26.9901.1033 (with 16" rear wheels) / 26.9901.3024 (with 20" to 24" rear wheels and hand rim)			

* Seat height adjustable by 10 cm depending on wheel size / ** Depending on wheel size



You can find more information on accessories in the Price List.

