

GMFCS
Levels IV-V



I am your universal seat pan sub-frame.

I carry you in almost all conventional seat pans and accompany you to many of your favourite locations – whether at home, in a facility or even, if you convert me to become a combination sub-frame, outdoors. The great thing about me is that you can adjust my seat angle and even my seat height to your requirements within seconds. Simply press the foot pedal, and I will adjust to floor level or to the height of the table.



I can also be converted into a combination sub-frame so that I can accompany you outside on the school playground or in the park.



**Seat angle**

can be adjusted to active and passive position.

Push handle

can be adjusted in angle.

Seat angle release lever

can also be mounted onto the push handle or the back of the pan.

**Trapezoidal adapter**

for all conventional seat pans.

Seat

can be lowered to a height suitable for play.

Foot pedal

adjusts the height of the sub-frame.

**Castors**

for easy transfer - with locking devices at the rear to prevent inadvertent rolling away.

Alternative: using a foot pump,

the height can be adjusted hydraulically.

**Frame colours**

- Lazurite blue metallic
- Pure white
- White aluminium
- Anthracite metallic
- Dormant red



Data	Size 0	Size 1	Size 2
Height of seat pan support	18 - 55 cm	18 - 55 cm	18 - 55 cm
Height of seat pan support with electrical height adjustment	26 - 55 cm	26 - 55 cm	26 - 55 cm
Height of seat pan support with hydraulic height adjustment	20 - 56 cm	20 - 56 cm	20 - 56 cm
Height of seat pan support with combination sub-frame	-	25 - 64 cm	25 - 64 cm
Length of base frame	58 cm	61 cm	75 cm
Total length with combination sub-frame	71 cm	76 cm	86 cm
Total width	52 cm	54 cm	64 cm
Total width with combination sub-frame	62 cm	66 cm	76 cm
Frame height	27 cm	27 cm	27 cm
Seat angle	(-)5° - 35°	(-)5° - 35°	(-)5° - 35°
max. load	30 kg	50 kg	75 kg
Weight	8 kg	10 kg	12 kg
HMV No.	26.9901.0022		

* Seat height adjustable by 8 cm depending on wheel size / ** Depending on wheel size



You can find more information on accessories in the Price List.

