

Features

- Manufactured by ISO7176
- Ultra-light weight: N.W. 26.6KGs , makes portability and transportation easy and convenient.
- Supporting passenger weighting up to 115kgs
- Areospace alumium and well-contructed frames achieves the solidity and safety
- Powful and safe Lithium Battery, 12V x 2 - 10Ah
- Powderful , dedicated 24V brushless motor drive

We are Offering

- Unique , excellent quality products with multiple models
- Competitive price
- Professional customer services
- Good warranty terms

Company's Certificates



Patent, Innovative, Multifunction Enrich Your ♥ Life

Unique Ultra-light Weight
Power Wheelchair
Easy to Maneuver



Innovative & Unique Power Wheelchair



4001

5001

3004

1002

Suzhou KD Intelligent Device Co., Ltd
No. 36, Gugang Road, Chenxiang Town
Taicang City, Jiangsu Province 215400

Tel: 0512-53110088 Ext 8001
Fax: 0512-53110099
E-mail: yiwenz@yahoo.com

Ultra-light weight-**26.6**kgs
Support weighing up
to **115**kgs

Powerful-dedicated 24V
Motor Drive
Lithium Battery 12Vx2-10AH

Cruise mileage **15**KM

Easy fold to Cars & SUV



Making the Life Easy and Convenient

Innovative and Unique Power Wheelchair



Special designed packaging or Travel case that make it well protected during transportation



Foldable in seconds, the answer for portability



Easy to fold & extremely convenient for traveling on any type of vehicle whether it is a Car, Truck, SUV, Taxi, Bus, Train, Cruise Ship or Airplane.



Length	Width	Height	Seat Depth	Seat Width	Seat Height	Backrest Height	N.W.
1000mm	660mm	860mm	420mm	380mm	525mm	410mm	26.6kg