

PL001-4001

## Features

- Manufactured by ISO7176
- Ultra-light weight: N.W. 30.5KGs , makes portability and transportation easy and convenient.
- Supporting passenger weighting up to 150kgs
- Areospace alumium and well-constructed frames achieves the solidity and safety
- Powful and safe Lithium Battery, 12V x 2 - 16AH
- Powderful , dedicated 24V brushless motor drive

## We are Offering

- Unique , excellent quality products with multiple models
- Competitive price
- Professional customer services
- Good warranty terms

## Company's Certificates



# Patent, Innovative, Multifunction Enrich Your ♥ Life

Unique Ultra-light Weight  
Power Wheelchair  
Easy to Maneuver



## Innovative & Unique Power Wheelchair



4001

5001

3004

1002

Suzhou KD Intelligent Device Co., Ltd  
No. 36, Gugang Road, Chenxiang Town  
Taicang City, Jiangsu Province 215400

Tel: 0512-53110088 Ext 8001  
Fax: 0512-53110099  
E-mail: yiwenz@yahoo.com

Ultra-light weight-**30.5**kgs  
Support weighing up  
to **150**kgs

Powerful-dedicated 24V  
Motor Drive  
Lithium Battery 12Vx2-16AH

Cruise mileage **18**KM

Easy fold to Cars & SUV



# Making the Life Easy and Convenient

Innovative and Unique Power Wheelchair



Mounted on Left or Right Side Joystick

Comfortable armrest

4 Position Adjustable Backrest

First Class safe Lithium Iron Phosphate battery, 12V x 2-16AH

Adjustable Leg Guard

Removable Cushion

Powerful Dedicated Brushless Motor Drive, 24v x 180w  
10" Rear Rubber Wheel

Storage Bag

Anti-tipping Wheels

Adjustable Footrest

Heavy duty 7" Front Wheel with Suspension & High End CNC Front Fork



Foldable in seconds, the answer for portability



Folded Size only **340x700x770mm**

Easy to fold & extremely convenient for traveling on any type of vehicle whether it is a Car, Truck, SUV, Taxi, Bus, Train, Cruise Ship or Airplane.



Length	Width	Height	Seat Depth	Seat Width	Seat Height	Backrest Height	N.W.
1100mm	700mm	920mm	420mm	420mm	510mm	450mm	30.5kg