



# HS-1500

*Connectorless Design, Combined Battery and Charger, Extra Space*



Outshining most normal power chairs, HS-1500 combines traditional functions and new ideals; adopting the Easy Connectorless Design, HS-1500 offers great convenience in maintenance. Users are now trouble free when dismantling the power chair and there is no complicated assembling procedure. And, with the combined battery and charger, charging HS-1500 is now as easy as one-two-three.

**CHIEN TI ENTERPRISE CO., LTD.**





## Features

- Maximum Speed : 6.5 km / 4 mph
- Capacity : 48.1 kg / 106 lbs
- Lightweight
- Durable Structure
- Connectorless Battery
- Height Adjustable /Liftable Armrest
- Left / Right Side Exchangeable Controller
- Length Adjustable Controller Assembly.

## TECHNICAL SPECIFICATIONS

Model Number	HS-1500
Overall Length	810 mm / 31.9"
Overall Width	590 mm / 23.2"
Overall Height	935 mm / 36.8"
Wheels: Front	130 mm / 5"
Wheels: Rear	200 mm / 8"
Weight w/ Batteries	48.1 kg / 106 lbs
Max. Speed	6.5 kmph / 4 mph
Ground Clearance	77 mm / 3"
Grade Climbable	6"
Curb Climbing	71 mm / 2.8"
Turning Radius	693.42 mm / 27.3"
Brake	Electro-Mechanical
Seat Type	Light Weight Foldable
Seat Width	460 mm / 18"
Drive Train	2-Motor, Rear-Wheel Drive
Motor Size	3600 rpm, 180 W
Battery Pack Weight	10.6 kg / 23.4 lb total
Travel Range	16 km / 10 mile
Battery Charger	2A
Electronics	PG VSI

\*Subject to change without notice.



## COLOR CHOICES



Blue



Burgundy



Silver



*Excellence in Mobility*

## CHIEN TI ENTERPRISE CO., LTD.

No.13, Lane 227, Fu Ying Rd., Hsin Chuang, Taipei, Taiwan,

Tel : 886-2-2903-2987

Fax : 886-2-2904-1170 · 2903-8807

E-Mail:sales@chienti.com.tw

<http://www.chienti.com.tw>